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Venduto a \_\_\_\_\_

Via \_\_\_\_\_ Località \_\_\_\_\_

Hangzhou Sejoy Electronics & Instruments Co. Ltd Building 2, No.202, Zhenzhong Road,  
West Lake Economy & Technology Zone 310030 Hangzhou China

EC REP

Shanghai International Holding Corp. GmbH (Hamburg) - Eiffelstrasse 80, 20537 Hamburg, Germany

Distrib. MORETTI S.P.A. Via Bruxelles, 3 - Melegnano 52022 Cavriglia (Ar) Italy - [www.morettispa.com](http://www.morettispa.com)

## INDEX

1. CODE .....	pag.
2. INTRODUCTION .....	pag.
3. INTENDED USE.....	pag.
4. DECLARATION OF CONFORMITY .....	pag.
5. GENERAL WARNINGS .....	pag.
6. ELECTROMAGNETIC COMPATIBILITY WARNINGS.....	pag.
7. SIMBOLOGY .....	pag.
8. GENERAL DESCRIPTION .....	pag.
8.1 Parts view .....	pag.
8.2 LCD Display .....	pag.
8.3 Basic Functions .....	pag.
9. USE WARNINGS .....	pag.
10. ADVANTAGES.....	pag.
10.1 Body temperature .....	pag.
11. HOW TO USE.....	pag.
11.1 Real Time Clock Setting .....	pag.
11.2 Get the room temperature.....	pag.
11.3 To measure ear temperature .....	pag.
11.4 To measure object temperature .....	pag.
11.5 After measurement .....	pag.
11.6 Temperature Taking Hints.....	pag.
11.7 Memory mode .....	pag.
11.8 Selecting a Temperature Scale.....	pag.
11.9 Battery Replacement.....	pag.
12. MAINTENANCE .....	pag.
12.1 Calibration.....	pag.
13. CLEANING AND DISINFECTION .....	pag.
13.1 Cleaning.....	pag.
13.2 Disinfection.....	pag.
14. CONDITIONS OF DISPOSAL .....	pag.
15. SPARE PARTS AND ACCESSORIES .....	pag.
16. TROUBLESHOOTING .....	pag.
17. ELECTROMAGNETIC COMPATIBILITY DECLARATION .....	pag.
18. TECHNICAL FEATURES.....	pag.
19. WARRANTY .....	pag.
20. REPAIRING .....	pag.
21. SPARE PARTS.....	pag.
22. EXEMPT CLAUSES.....	pag.

IIa Class Medical Device

According to Directive 93/42/EEC and further modifications

## 1. CODE

RM355 - Infrared ear thermometer

## 2. INTRODUCTION

Thank you for purchasing a GRADO a by Moretti Thermometer. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance. In case of doubts, we recommend to contact the dealer who can help and advice you properly.

### CAUTION!

- Do not use the product for a purpose not indicated in this manual.
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- The manufacturer reserves the right to change the information contained in this document without previous notice.

## 4. DECLARATION OF CONFORMITY

Moretti S.p.A. declares under his sole responsibility that the product placed in the market by the same Moretti S.p.A. is in conformity to the relevant provisions of Directive 93/42/CEE on medical devices and subsequent amendments and are produced by Hangzhou Sejoy Electronics & Instruments Co. Ltd according to the standards indicated in the regulations:

- CEI 60601-1 - Part 1: General requirements for basic safety and essential performance
- EN 60601-1-2 - Part 1: General requirements for safety - Collateral standard: Electromagnetic compatibility
- ISO 80601-2-56:2009:Medical electrical device - (Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for measuring body temperature).
- The manufacturer is ISO 9001 and ISO 13485 certified.

## 5. GENERAL WARNINGS

- Read instructions thoroughly prior to use.
- Store and use thermometer away from mobile phones.
- Always take temperature readings from the same ear, as ears may differ in temperature.
- Keep the battery away from children.
- Keep the unit dry and away from any liquids and direct sunlight.
- Infrared thermometers should not be used shortly after exercise, bathing or coming indoors.
- Choking Hazard: Thermometer cap and battery may be fatal if swallowed.  
Do not allow children to use this device without parental supervision.
- Do not store or keep in direct sunlight, or in high temperature (below -25°C/-13°F or over 55°C /131°F) non excessive humidity (>95%RH). The packaging is made in paper.
- Shelf life - the thermometer use limit is defined by the wear of the parts.
- Do not place thermometer battery near extreme heat as it may explode.
- Do not allow children to play on or operate the thermometer.
- It is recommended the performance should be checked every two years.
- Remove battery from the device when not in operation for a long time.
- The use of temperature readings for self-diagnosis may lead to the worsening of existing disease conditions.
- Do not attempt measurements when the thermometer is wet as inaccurate readings may result.
- Do not bite the thermometer. Doing so may lead to breakage and/or injury.
- Do not attempt to disassemble or repair the thermometer. Doing so may result in inaccurate readings.

- After each use, disinfect the thermometer especially in case the device is used by more than one person.
- This thermometer must only be used with this brand thermometer probe covers.
- To avoid inaccurate readings, always attach a new, clean probe cover for each temperature reading.

## 6. ELECTROMAGNETIC COMPATIBILITY WARNINGS

This equipment has been tested and it complies with medical devices limits belonging to EN 60601-1-2:2007. These limits have been designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guaranteed that the interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference following one or more of the suggestions below:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Consult the manufacturer or field service technician for help.
- Do not use in proximity of the device, mobile phones and other devices, which can generate strong electrical or electromagnetic fields as they can interfere with the operation of the device.

## 7. SIMBOLOGY


 serial number  warning  EC REP  type BF  instruction

 lot  date of production  disposal  CE Mark  manufacturer


## 8. GENERAL DESCRIPTION

The Infrared Ear thermometer is specifically designed for safe use in the eardrum. The Infrared Ear Thermometer is a device capable of measuring people's body temperature by detecting the intensity of infrared light emitted from the ear canal of human. It converts the measured heat into a temperature reading displayed on the LCD. When properly used, it will quickly assess your temperature in an accurate manner.

### 8.1 Parts view

## 8.2 LCD Display

## 8.3 Basic Functions

### Real Time Clock

The real time clock will be recorded with the memory function and help you to recognize each measurement result.

- Please see the Real time clock setting section to learn how to setup the time in the first use.

### Room Temperature

Suitable ambient temperature is important for the baby and patient, the thermometer always help you recognize the room temperature.

- Please see the Illustration For Use section to learn how to get the correct room temperature.

### Ear Temperature

Mode The thermometer has been designed for practical use. It's not meant to replace a visit to the doctor. Please also remember to compare the measurement result to your regular body temperature.

- Please see the Illustration For Use section to learn how to measure the body temperature.

### Object Temperature

Mode The object mode shows the actual, unadjusted surface temperatures, which is different from the body temperature. It can help you to monitor if the object temperature is suitable for the baby or patient, for example the baby's milk.

- Please see the Illustration For Use section to learn how to measure the object temperature.

### Fever Alarm

If thermometer detects a body temperature  $\geq 37.8^{\circ}\text{C}$  under ear mode, there will be a long beep sound followed by three short beep sound to warn the user for potential fever.

### Memory Mode

There are each 10 sets memories for ear and object measurements. Each memory also records the measurement date/time/mode icon.

### C / F Switch

Please see the Selecting A Temperature Scale section to learn how to change between Celsius and Fahrenheit.

### Probe Cover

Automatic Detection Automatic detection function, convenient for test with or without probe cover.

## 9. USE WARNINGS


The device performance can be compromised if one or more of the following events occur:

- Operation outside the manufacturer's stated temperature and humidity range.
- Storage outside the manufacturer's stated temperature and humidity range.
- Patient temperature is below ambient temperature
- Do not bend, drop the thermometer.
- Do not boil probe to disinfect.
- Do not disassemble, except in case of battery replacement.
- The thermometer should only be used under the supervision of an adult. Prohibit children from biting, hitting and removing the battery cover, dislodging and potentially swallowing the battery.
- Do not use the device in case of strong electromagnetic field such as microwave oven or cell phone.

## 10. ADVANTAGES

Infrared Ear Thermometer measures core body temperature, which is the temperature of a body's vital organs.

### Vital Organs

(See Figure 3). Ear temperatures accurately reflect core body temperatures since the eardrum is in the cranial cavity and is not affected by the environment. The eardrum shares blood supply with the hypothalamus (part of the brain.) Therefore, body temperature changes are reflected sooner in the ear than they are in other parts of the body.

Clinical studies have shown that the ear is an excellent site for temperature measurement because temperatures taken in the ear reflect the body's core temperature. Body temperature is regulated by the hypothalamus, which shares the same blood supply as the tympanic membrane. Changes in core body temperature are usually seen sooner at the tympanic membrane than at other sites, such as the rectum, mouth or under the arm.

Advantages of taking temperatures at the ear versus other sites:

**Fig.3** **Axillary:** temperature readings only reflect skin temperature which may not indicate the internal body temperature.

**Rectal:** temperatures often lag significantly behind internal body temperature changes, especially at times of rapidly changing temperatures.

**Oral:** temperatures are often influenced by eating, drinking, thermometer placement, breathing through the mouth, or the inability of the person to close their mouth completely.

### 10.1 Body temperature

- Body temperature can vary from 35°C to 38°C. The generally accepted "normal" temperature is 37.0°C.
- Usually your temperature is lower in the morning than in the afternoon. Normal temperature varies according to the location of the temperature reading.

Normal body temperature is a range. The following table shows that ranges of normal also vary by site. Therefore, readings from different body locations, even taken at the same time, should not be directly compared.

Normal ranges by body locations:

Axillary	34.7 - 37.3°C	94.5 - 99.1°F
Oral	35.5 - 37.5°C	95.9 - 99.5°F
Rectal	36.6 - 38.0°C	97.9 -100.4°F
Ear	35.8 - 38.0°C	96.4 -100.4°F

Also, a person's normal temperature range tends to decrease with age.


The following table shows normal ranges by age:

0-2 Years	36.4 - 38.0°C	97.5 - 100.4°F
3-10 Years	36.1 - 37.8°C	97.0 - 100.0°F
11-65 Years	35.9 - 37.6°C	96.6 - 99.7°F
Over 65 Years	35.8 - 37.5°C	96.4 - 99.5°F

The range of normal varies from person to person and can be influenced by many factors such as time of day, level of activity, medications and gender.

## 11. HOW TO USE

### 11.1 Real Time Clock S etting

When use thermometer for the first time, please set the parameters of the thermometer. With the thermometer off, press and hold  BUTTON to enter into setting mode for 3 seconds.

#### 1. Set the time format

The device can display the time in either an AM/PM (12-hour) or a 24:00 (24-hour) format. Press and release TEST BUTTON to select the format. With the preferred time format on the display, press  BUTTON, the Hour figure is flashing automatically.

#### 2. Set the hour

Press and release the TEST BUTTON to advance one hour until the correct hour appears. After the hour is set, press  BUTTON, the Minute figure is flashing automatically.


#### 3. Set the minute

Press and release the TEST BUTTON to advance one minute until the correct minute appears. After the minute is set, press  BUTTON, the Year figure is flashing automatically.

#### 4. Set the year

Press and release the TEST BUTTON to advance one year until the correct year appears. After the year is set, press  BUTTON, the Month figure will appear.

#### 5. Set the month


Press and release the TEST BUTTON to advance one month until the correct month appears. After the month is set, press  BUTTON, the Date figure is flashing automatically.

#### 6. Set the date

Press and release the TEST BUTTON to advance one day until the correct date appears. After the date is set, press  BUTTON to exit the setting mode.

## 11.2 Get the room temperature

1. When the thermometer has been turned on, room temperature will be shown on the screen with icon. Room temperature, date and time will alternately display every 3 seconds. (See figure 4)
2. For room temperature, the thermometer should be placed on the table and avoid the direct sunshine or air conditioner flow on.
3. The room temperature can be reference after 5 minute later when taking the temperature.

(Fig.5)

## 11.3 To measure ear temperature

1. To achieve accurate readings, make sure a new, clean probe cover is in place before each measurement. Place a new probe cover on the probe firmly; making sure the probe and the inner portion of the probe cover fit snugly. To ensure proper hygiene and to avoid contamination, a new probe cover should be applied for each user.
2. Press the TEST BUTTON, you can get the body temperature by ear measurement after hearing two beep sounds.

(Fig.6)

## Remarks:

- a. Probe cover sign stand for taking the ear temperature with probe cover now and vice versa.
  - b. To avoid the risk of cross contamination, it is recommended that you always measure the ear temperature with probe cover.
3. If you see as shown in figure 6 shows, it stand for the mode of taking the object temperature. You can press and hold  BUTTON and then press and hold the TEST BUTTON to switch the mode between ear mode and object mode.
  4. Insert the probe along the ear canal as deep as possible into the ear as depicted in Figure 7. Then press the TEST BUTTON once to start the measurement. A long beep will sound approximately 1 second after the TEST BUTTON is activated. This signals the end of the process as results are displayed on the LCD screen.
  5. The thermometer is ready for testing again once an ear mark displays on the LCD. A waiting period of 20 seconds between testing is recommended to avoid excessive cooling of the skin.
  6. The thermometer will shut off automatically after 60 seconds of inactivity. To prolong battery life, press the  BUTTON to turn the unit off.

## 11.4 To measure object temperature

(Fig.8)

1. Eject the probe cover out firstly. Press the TEST BUTTON to turn on the thermometer, you can take the object temperature after hearing two beep sounds. (see figure 5)

## Remarks:

- a. This mode shows the actual, unadjusted surface temperatures, which is different from the body temperature.
- b. If you see as shown in figure 3 shows, it stand for the mode of taking the ear temperature. You can press and hold  BUTTON and then press and hold the TEST BUTTON to switch the mode between ear mode and object mode.

c. If "Er6" occurs, please refer to the section of for details.

### 11.5 After measurement

1. Power off: Device will automatically shut off if left idle for more than 1 minute to extend battery life.
2. Clean the probe after each use to ensure an accurate reading and avoid cross contamination. (See the section of Care and Cleaning for details.)

### 11.6 Temperature Taking Hints

1. The right ear reading may differ from the reading taken at the left ear. Therefore, always take the temperature in the same ear.
2. The ear must be free from obstructions or excess earwax buildup to take an accurate reading.
3. External factors may influence ear temperatures such as:
  - been lying on one ear or the other
  - had their ears covered
  - been exposed to very hot or very cold temperatures
  - been recently swimming or bathing

In these cases, remove the individual from the situation and wait 20 minutes prior to taking a temperature.

4. For persons wearing hearing aids or ear plugs, remove the device and wait 20 minutes prior to taking a temperature.

(Fig.10)

(Fig.9)

5. When using the thermometer on infants under age 1, pull the ear up making sure the sensor faces the eardrum. (See Figure 9)
6. When using the thermometer on individuals over the age of 1, pull the ear back making sure the sensor faces the eardrum. (See Figure 10)
7. Do not touch the probe window. A cotton swab and rubbing alcohol may be used to clean the sensor window.
8. If thermometer is stored in a significantly different environment than testing location, place it in the testing location for approximately 30 minutes prior to use.

### 11.7 Memory Mode

(Fig. 11)

(Fig. 12)

1. The Memory Mode can be accessed either in ear mode or object mode:
  - When the thermometer has been turned on and followed by Figure 3 or finished testing, press and hold the  Button for three seconds. The letter M will appear in the center left corner of the display. (See Figure 8)
2. The thermometer will automatically memorize the last 10 temperature readings. Each memory also records the measurement date/time/mode icon. Each time the  Button is pressed, the screen displays past readings that correspond with a number 1-10. The number 1 reflects the most recent reading while the number 10 reveals the oldest reading stored in memory. (See Figure 9)
3. In the memory mode, ear mark or scan mark always exist. The user can press the TEST BUTTON to take new measurements.


### 11.8 Selecting A Temperature Scale

1. Temperature readings are available in the Celsius (°C) or Fahrenheit (°F) scale.
2. With the unit off, press and hold the TEST BUTTON for 3 seconds to enter into unit changing mode.
3. Press and release TEST BUTTON to select the unit.
4. When the preferred unit on the display, press  BUTTON to exit the unit changing mode.

### 11.9 Battery Replacement

1. Replace battery when " " appears in the upper right corner of LCD display. (See Figure 11)

- Slide battery cover down as shown in Figure 12.
- Remove battery and install 2 new AAA alkaline batteries as shown in Figure 13 .
- Slide battery cover back on.

(Fig. 11)

## 12. MAINTENANCE

The GRADO devices by Moretti S.p.A are checked carefully and supplied with EC mark, once launched on the market. For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every 2 years. In case of repair, use only original spare parts and accessories.

### 12.1 Calibration

The thermometer is initially calibrated at the time of manufacture. If the thermometer is used according to the use instruction, periodic readjustment is not required. However, We recommends checking calibration every two years or whenever clinical accuracy of the thermometer is in question. Please send the complete device to the dealers or manufacturer. The above recommendations do not supersede the legal requirements. The user must always comply with legal requirements for the control of the measurement, functionality, and accuracy of the device which are required by the scope of relevant laws, directives or ordinances where the device is used.

## 13. CLEANING AND DISINFECTION

### 13.1 Cleaning

- Probe covers are intended for single use only!  
Absent, damaged, perforated, or soiled probe covers can lead to inaccurate temperature readings.
- The probe window must be kept clean, dry, and undamaged at all times to ensure accurate readings. The accuracy of temperature readings can be affected by damage to the probe window, or the presence of dirt and ear wax on the probe window.
- Fingerprints, earwax, dust and other soiling compounds reduce transparency of the window and result in lower temperature readings. If the thermometer is accidentally placed in the ear without a probe cover installed, the window must be cleaned immediately.
- The probe window is the most delicate part of the thermometer. To safely clean the window, gently wipe its surface with a cotton swab slightly moistened with isopropyl alcohol and immediately wipe dry with a clean cotton swab. After cleaning, allow at least 5 minutes drying time before taking temperatures.

**Note:** Do not use any chemical other than isopropyl alcohol to clean the probe window.

- Use a soft, dry cloth to clean the thermometer display and exterior.
- The thermometer is not waterproof. Do not submerge the unit in water when cleaning.
- Store thermometer and probe cover in a dry location, free from dust and contamination and away from direct sunlight.
- Periodic cleaning and disinfection of the device following use to prevent patient cross infection.
- Use a soft cloth slightly moistened with a 70% isopropyl alcohol solution to disinfect the thermometer and probe. Do not use abrasive cleaners.
- Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.
- Do not remove or modify the equipment without permission.
- Strong electromagnetic fields may interfere with the proper operation of the thermometer. Do not use this thermometer close to sources of strong electromagnetic radiation.

### 13.2 Disinfection

If you need to disinfect the device use isopropyl rubbing alcohol

**NOTE:** Never use acids, alkalis or solvents such as acetone

## 14. CONDITIONS OF DISPOSAL

### General conditions of disposal

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the thermometer in the appropriate disposal areas for recycling.

WARNINGS FOR THE PROPER DISPOSAL OF THE PRODUCT, AS REQUIRED BY THE EUROPEAN DIRECTIVE 2002/96/EC: (DIRECTIVE 2002/96/EC)

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin the product bears the barred bin symbol.

DISPOSAL OF WASTE BATTERIES (DIRECTIVE 2006/66/EC)

This symbol on the battery or on the packaging indicates that the battery provided with this product shall not be treated as household waste. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. The recycling of the materials will help to conserve natural resources. At the end of their life hand the batteries over to the applicable collection points for the recycling of waste batteries. For more detailed information about recycling of this product or battery, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

## 15. SPARE PARTS AND ACCESSORIES

For spare parts and accessories refer to main catalogue.

RM360 - Probe cover

## 16. TROUBLESHOOTING

Problem	Solution	Error message
Measurement before thermometer is ready	Take a measurement until ORECCHIO or WI_FI appears on the display	Er 1
The ambient temperature is not within the range between 10°C and 40°C (50°F~104°F).	Place the thermometer in a room for at least 30 minutes at room temperature between 10°C and 40°C (50°F~104°F)	Er 2
The thermometer is placed incorrectly or unsteady	Read Illustration For Use thoroughly and take a new temperature measurement	Er 3
The thermometer showing a rapid ambient temperature change	Allow the thermometer to rest in a room for at least 30 minutes at room temperature: between 10°C and 40°C (50°F~104°F)	Er 4
The thermometer is not functioning properly.	Unload the battery, wait for 1 minute and repower it. If the message reappears, contact the retailer for service.	Er 5

In object temperature mode, measurement with probe cover.	Remove probe cover in object temperature mode.	
In Ear mode: Temperature taken is higher than 43.0 °C (109.4°F). In Object mode: Temperature taken is higher than 100°C (212°F).	Read Temperature Taking Hints Thoroughly, then check the integrity of the probe cover and take a new temperature measurement.	
In Ear mode: Temperature taken is lower than 28.0°C (82.4°F). In Object mode: Temperature taken is lower than -20°C (-4°F). Read Temperature	Taking Hints thoroughly , then make sure the probe cover and lens filter are clean, then take a new temperature measurement.	
The thermometer could not work due to low battery.	Replace two new alkaline batteries size AAA	

### 17. ELECTROMAGNETIC COMPATIBILITY DECLARATION

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The device uses RF Energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR11	Class B	The device is suitable for use in all establishment. Including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	N/A	
Voltage fluctuations/flicker emissions IEC61000-3-3	N/A	

### Guidance and manufacturer's declaration - Electromagnetic immunity -

The devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	± 2kV for power supply line ± 1kV for input/out line	N/A	
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	N/A	

Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % UT (>95 % dip in UT)for 0,5 cycle 40 % UT (60 % dip in UT)for 5 cycles 70 % UT (30 % dip in UT)for 25 cycles <5 % UT (>95 % dip in UT)for 5 sec	N/A	
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment.
<b>Note: UT is the a.c. mains voltage prior to the application of the test level</b>			

#### Guidance and Manufacturer's Declaration - Electromagnetic Immunity -

The devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3Vrms150 kHz to 80 MHz outside ISM bandsa	N/A	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 80 MHz to 2.5G MHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey c, should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

NOTE 3 Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and landmobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device

NOTE 4 Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 KHz to a 80MHz $d=1.2\sqrt{P}$	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2,5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73

1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

The user can prevent electromagnetic interference by maintaining a minimum distance between the thermometer and the device that might causes interference.

## 18. TECHNICAL FEATURES

<b>Measuring site:</b>	Ear canal, Room, objects, liquids
<b>Measuring range:</b>	Ear Mode: 28.0°C ~ 43.0°C (82.4°F ~ 109.4°F) Object mode: -20°C ~ 100°C (-4°F ~ 212°F)
<b>Room temperature:</b>	0°C ~ 50°C (32°F ~ 122°F)
<b>Laboratory accuracy:</b>	Ear mode: +0.2°C (0.4°F) during 35.5°C ~ 42.0°C (95.9°F ~ 107.6°F) at 15°C ~ 35°C (59.0°F ~ 95.0°F) operating temperature range +0.3°C (0.5°F) for other measuring and operating temperature range Object mode: +4% or +2°C (4°F) whichever is greater Room temperature: +2°C (+4°F)
<b>Display:</b>	Liquid crystal display (LCD) with backlight 0.1°C or 0.1°F
<b>Response time:</b>	1 seconds
<b>OPERATING</b>	
<b>Temperature range:</b>	10.0°C ~ 40.0°C (50°F ~ 104°F); 15% ~ 85% RH, non-Condensing
<b>Storage and transportation</b>	
<b>Temperature range:</b>	-25°C to 55°C (-13°F to 131°F); 15% ~ 95% RH, non-Condensing
<b>Clinical accuracy:</b>	Clinical bias: 0.04°C (0.07°F) Clinical repeatability: 0.11°C (0.2°F) Limits of agreement: 0.69°C (1.2°F) +0.3°C (0.5°F)
<b>Dimension (lxwxh):</b>	155x38x47 mm
<b>Weight</b>	
<b>(Including battery):</b>	Approx 83gr (with batteries)
<b>Battery:</b>	3V DC (2xAAA Batteries)
<b>Battery life:</b>	Approx 1 Year/3000 readings
<b>Classification:</b>	tipo BF 
<b>Ingress protecting</b>	
<b>Rating:</b>	IP22

## 19. WARRANTY

Every Moretti product is guaranteed for 2 (two) YEARS from the purchase date, except the restrictions specified below. The warranty shall not apply in case of damage caused by improper use, abuse, alterations of the product or any other use not in compliance with the user manual. The proper use of the product is indicated in the user manual. Moretti S.p.A. shall not be liable for any damages or personal injuries due to incorrect assembly/improper use of the device, not compliant with the user manual. Moretti does not guarantee its products against damages or faults in the following cases: natural disasters, non authorized maintenance or reparations, damages caused by electrical supply problems (if provided), use of parts or components not provided by Moretti, failure in accomplishing the guidelines and instructions for use, non authorized alterations, damages due to the transport (different from Moretti original delivery) or failure in maintaining the device as indicated in the manual. The components subject to wear are not covered by this warranty if the damage is due to the normal use of the product.

## 20. REPAIRING

### Warranty repair

If a Moretti product shows material or manufacturing faults during the warranty period, Moretti S.p.A., together

with the client, will consider whether the product fault is covered by warranty. Moretti S.p.A., at its unquestionable discretion, may replace or repair the item under warranty in a Moretti reseller point or at its own qualified offices. Labour costs for repairing the product may be at Moretti's expense if it is determined the reparation is under warranty. A reparation or replacement does not renew nor extend the warranty.

#### Repair not covered by warranty

A product, not covered by warranty, may be returned for reparation only after receiving prior authorization from Moretti client service. Labour and delivery costs for an out-of-warranty reparation will be fully at the client's or reseller's expense. Reparations of out-of-warranty products are guaranteed for 6 (six) months, from the day the client receives the repaired product.

#### Non defective device

The client is informed if, after examining and testing a returning product, Moretti decides the product is not defective. The product will be sent back to the client and all the return delivery costs will be at his expense.

## 21. SPARE PARTS

The original Moretti spare parts are guaranteed for 6 (six) months from the day you receive them.

## 22. EXEMPT CLAUSES

Except what clearly specified in this warranty and in accordance with the law, Moretti S.p.A. does not offer its clients any other declaration, warranty or condition, expressed or implicit, including any declaration, warranty or condition of merchantability, of non violation or non interference, of suitability for a specific aim. Moretti S.p.A. does not guarantee that the use of Moretti products will be continuous and faultless. The duration of possible implicit warranties that may be applied in accordance with the law is limited to the warranty period.

Some states or countries do not allow limitations to the duration of an implicit warranty, or the exoneration, or the limitations of accidental or indirect damages related to products for consumers. In such states or countries some exoneration or limitation of this warranty may not be applied to the user. This warranty is subject to modifications without any prior notice.

### WARRANTY CERTIFICATE

Product \_\_\_\_\_

Purchased on (date) \_\_\_\_\_

Retailer \_\_\_\_\_

Address \_\_\_\_\_ Town/city \_\_\_\_\_

Sold to \_\_\_\_\_

Address \_\_\_\_\_ Town/city \_\_\_\_\_

Hangzhou Sejoy Electronics & Instruments Co. Ltd Building 2, No.202, Zhenzhong Road,  
West Lake Economy & Technology Zone 310030 Hangzhou China

Shanghai International Holding Corp. GmbH (Hamburg) - Eiffestrasse 80, 20537 Hamburg, Germany

Distrib. MORETTI S.P.A. Via Bruxelles, 3 - Meleto 52022 Cavriglia (Ar) Italy - [www.morettispa.com](http://www.morettispa.com)

**MORETTI S.P.A.**

Via Bruxelles, 3 - Meleto  
52022 Cavriglia (Arezzo)

Tel. +39 055 96 21 11  
Fax. +39 055 96 21 200

[www.morettispa.com](http://www.morettispa.com)  
[info@morettispa.com](mailto:info@morettispa.com)